

# Birdie 150 Compact

COMPACT AND COMFORTABLE



Designed to fold easily and without the need for tools, the Birdie Compact can be stored or transported comfortably when space is limited. The compact base and large castors manoeuvre smoothly around beds, in bathroom areas and through narrow doorways.



Yes, you can.®

Making Life's  
Experiences Possible

# Birdie 150 Compact



The Birdie Compact is only 55 cm wide making it the only lifter that is TÜV certified that is this narrow - unique to the market. The ergonomically designed push bar allows the carer to position their hands comfortably reducing the likelihood of strain. The 100mm smoothing rolling castors make it easy to push and pull the lifter.

The Birdie Compact takes up minimal space when stored reducing clutter in confined spaces. When space is limited, the Birdie Compact can even be stored upright.

Users can enjoy maximum flexibility with a range of spreader bars and a very comfortable sling range, which allows for individual adjustment of the transfer. The quick release attachment system makes configuration even easier.

## Features



Quick release attachments.



Tool-less and easy to fold.



Intelligent monitoring system.



Store upright to save space.


## Tech Specs



**Length**  
109cm



**Min. base width**  
52cm



**Max. base width**  
84.5cm



**Arm height**  
52.5 - 161.5cm



**Ground clearance**  
10cm



**Max. user Weight**  
150kg



**Unit weight**  
36kg



**Battery time**  
40 cycles

### Warranty:

3 years frame  
2 years electronics and motors  
1 year battery



FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [sales@invacare.com.au](mailto:sales@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

**Asia:** Phone +61 2 8839 5330, email [exportsales@invacare.com.au](mailto:exportsales@invacare.com.au) or visit [www.invacare.com](http://www.invacare.com)

**Thailand:** Phone +66 0 2108 8555 or visit [www.invacare.com](http://www.invacare.com)



Yes, you can.®