

The ergonomically designed Pegasus METRO has been designed from the ground up to focus on key aspects important to both new and seasoned scooter users. As a result, it is equipped with numerous features to ensure users experience a safe, easy and enjoyable ride.





Pegasus Metro

The suspension system on the compact Pegasus METRO has been designed to navigate a variety of everyday surfaces and terrains, allowing you to relax and enjoy a controlled, sturdy drive in comfort.

In addition to a unique lighting system, the Pegasus METRO is equipped with an electronic speed reduction feature that automatically reduces speed around bends, ensuring safe, comfortable manoeuvres at all times. Furthermore, a two-step disengaging lever prevents the scooter from free-wheeling if it is knocked out of drive mode.

Electronics and the motor are protected from corrosion and water which safeguard the vital workings of the scooter against dirt preventing any unnecessary malfunctions.

Features



Users can easily adjust the tiller with a lever to suit their needs.



Ergonomically designed steering system makes for easier, more



Ensures immediate braking if required.



Clear, bright brake lighting warns people that the scooter is slowing down.



Dashboard configuration.

Tech Specs



Backrest width 510 mm



Seat depth 470 mm



Floor to seat height Min. 440 mm Max. 510 mm



Backrest height 490mm



Driving unit width 660 mm



Total length 1270 mm



Empty weight 110 kg



Max user weight 136 ka



Battery capacity 2 x 12 V/50 Ah



Turning diameter



Max climbable obstacle height 80 mm



Ground clearance





Max driving distance 54 km



Max driving speed



Motor capacity

Colour Options



Alium Blue



Iconic Chrome Diamond Silver Marble White Hexagon White









Warranty Terms and Conditions apply. Contact Invacare for full details.

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

